

# Omega Laboratories Oral Fluid Drug Panels & Cutoff Levels

	Screening (Neat)	Confirmation (Neat)	10 Panel <i>S10P</i>	9 Panel <i>S9P</i>	6 Panel <i>S6P</i>	5 Panel <i>S5P</i>
<b>Amphetamines</b>	300 ng/ml		●	●	●	●
Amphetamine		50 ng/ml	●	●	●	●
<b>Methamphetamine / MDMA</b>	120 ng/ml		●	●	●	●
Methamphetamine		50 ng/ml	●	●	●	●
MDA		50 ng/ml	●	●	●	●
<b>Cocaine / Metabolite</b>	15 ng/ml		●	●	●	●
Cocaine		8 ng/ml	●	●	●	●
Benzoylcegonine		8 ng/ml	●	●	●	●
<b>Marijuana (THC)</b>	3 ng/ml		●		●	
THC		2 ng/ml	●		●	
<b>Opiates</b>	30 ng/ml		●	●	●	●
Codeine		Morphine 40 ng/ml	●	●	●	●
Hydrocodone		40 ng/ml	●	●	●	●
6-MAM (Heroin)		4 ng/ml	●	●	●	●
<b>Phencyclidine</b>	3 ng/ml		●	●	●	●
Phencyclidine (PCP)		10 ng/ml	●	●	●	●
<b>Oxycodone</b>	30 ng/ml		●	●		
Oxycodone		40 ng/ml	●	●		
<b>Methadone</b>	15 ng/ml		●	●		
Methadone		15 ng/ml	●	●		
<b>Benzodiazepines</b>	3 ng/ml		●	●		
Alprazolam		Nordiazepam 5 ng/ml	●	●		
Diazepam		5 ng/ml	●	●		
<b>Barbiturates</b>	60 ng/ml		●	●		
Secobarbital		Butalbital 60 ng/ml	●	●		
Pentobarbital		Phenobarbital 60 ng/ml	●	●		
Amobarbital		60 ng/ml	●	●		
<b>Separate Panels</b>			(Offered separate from 6 and 10 panels)			
<b>Cotinine (Nicotine Metabolite)</b>	30 ng/ml		(Offered separate from 6 and 10 panels)			
Cotinine		10 ng/ml	(Offered separate from 6 and 10 panels)			

**Neat vs. Dilute:** Screening and Confirmation listed as Neat in accordance with recent SAMHSA guidelines.

**Validity Testing:** All oral fluid panels include IgG human validity testing